# Characterizing Peptide-Antibody Interactions for Disease Detection

Jane Lin

2<sup>nd</sup> Year Biology major with a concentration in Bioengineering

Mentor: Serra Elliott

Faculty Advisor: Dr. Patrick S. Daugherty


Dept. of Chemical Engineering

EUREKA!, UCSB

08/29/13


## Pre-eclampsia (PE) is One of the Top 5 Leading Causes of Maternal Death

- $\approx 5 8\%$  pregnant woman (1)
- Advanced conditions:
  - HELLP syndrome
  - Seizures
- Treatment: Delivery of baby





# Potential Role of Immunological Component in Pathogenesis





- Poor Diagnostic Tool
  - Clinical test is lengthy and complicated
  - Normal-outcome pregnancies & Other diseases


Hypothesis: Additional PEspecific antibodies exist and could serve as effective biomarkers

## Antibody Detecting Peptides for Diagnostics


- Fluorescence signal = peptide-antibody binding
  - peptide expression and target concentration
- Accuracy of a peptide-antibody interaction to differentiate between PE and Normal-outcome








### Fluorescent Signal Increases with Peptide Density




**Induction Times (Peptide Density)** 

#### Peptide Density Also Affects Apparent Binding Affinity



## Mean Fluorescent Fold of PE Patients Indicates Higher Reactivity with Peptide




<sup>\*</sup>t-test,\*\*U-test p < 0.05

# Receiver Operating Characteristic (ROC) Curve Determines Overall Peptide Diagnostic Efficacy

$$\frac{\text{Correct PE}}{\text{Total # of PE}} = 0.65$$

$$1 - \frac{\text{Correct Normal}}{\text{Total # of Normal}} = 0.15$$




## Peptide Has Potential Use as a Diagnostic Tool Specific to PE

- Direct correlation between peptide density and fluorescence intensity
- Peptide density affects cellsurface binding affinity
- Area under ROC curve indicates fairly high diagnostic efficacy



#### Future Prospects for PE Patients



#### Acknowledgments

Serra Elliott Dr. Patrick Daughtery Daugherty Group





